

SAM
MEDICAL

SAM[®] Pelvic Sling II

**SLIDE.
CLICK.
SECURE.**

OPTIMAL FORCE STABILIZATION

SAM[®] Pelvic Sling is the first and only force-controlled circumferential pelvic belt. Featuring an optimized one-piece design, it's scientifically proven to safely and effectively reduce and stabilize open-book pelvic ring fractures.



FEATURED TECHNOLOGY

AUTOSTOP BUCKLE

The patented SAM[®] AUTOSTOP buckle acts to avoid over- or under-tightening of SAM[®] Pelvic Sling. Once optimal compressive force is reached, two prongs activate, clicking to confirm correct application. The sling is secured in place with its hook and loop system strap.



ENGINEERED TO PRESERVE LIFE

1 LOW-FRICTION POSTERIOR SLIDER

Easier transfers following severe trauma made possible by low-friction materials covering backside of the pelvic sling.

2 SAM[®] AUTOSTOP BUCKLE

Tactile and audible “click” provides clear feedback, confirming proper application and preventing over-tightening while maintaining optimal force.

3 PRECISION STRAP


Precisely-engineered holes lock with AUTOSTOP Buckle, allowing for quick application. Secure fastening with sturdy hook and loop strap.

4 TAPERED DESIGN

Narrow and tapered front facilitates urinary catheterization, interventional radiology, external fixation, and abdominal surgery without removal.

5 RADIOLUCENT MATERIAL

All materials are radiolucent except the buckle’s two stainless steel springs, allowing X-rays and CT-scans without removal.



MODELS



CIVILIAN
ORANGE/BLUE



MILITARY
OD GREEN /BLACK

PRODUCT SPECIFICATIONS

PACKAGED SIZE & WEIGHT*

Height: 5.5 in (14 cm)
Width: 7.5 in (19 cm)
Depth: 3.75 in (9.5 cm)
Weight: 8.6 - 9.2 oz (244 - 250g)**

OPERATING PRINCIPLE

Optimal force is reached when 150N of circumferential pressure is attained, which extends prongs from SAM[®] AUTOSTOP Buckle through the Precision Strap.

MODELS AVAILABLE

Military
Standard 32-50"/ 81-127 cm (SL556652-OD)

Civilian
Small 27-45"/ 69-114 cm (SL556652-SM)
Standard 32-50"/ 81-127 cm (SL556652)
Large 36-54"/ 91-137 cm (SL556652-LG)

GLOBAL SALES

U.S. Toll Free: 800.818.4726
U.S. Local: 503.639.5474
U.S. Fax: 503.639.5425
sammedical.com



Ergo Europe, Prinsessegracht 20,
2514 AP, The Hague, The Netherlands



SAM Medical Products[®]
27350 SW 95th Ave, Ste 3038
Wilsonville, OR 97070 USA

* Size and weight vary by configuration and manufacturing process.

**Weight varies by size (Small, Standard, Large)

SAM[®] is a registered trademark in the United States and other countries.

©2017 SAM Medical Products. All rights reserved.