

Blist-O-Ban®

It takes a blister to prevent a blister. That's the science behind Blist-O-Ban®. Its patented BursaTek® dome glides smoothly in all directions, deflecting friction and shear forces away from the skin. Place the ultra-thin Blist-O-Ban over any hotspots, blisters or injury prone areas.

A study funded by the National Institute of Health demonstrated this new technology to be significantly more effective at reducing friction than ten other commonly used products. Temple University School of Podiatry performed a study involving 100 subjects who walked 20 miles in warm rainy conditions. These walkers developed 5 times more blisters in areas not protected by the Blist-O-Ban. Field-tests involving special operations soldiers and competitive athletes also confirmed these laboratory studies.



FEATURES

- Medical-grade adhesive
- Water-resistant
- Multi-day protection
- Latex-free
- Breathable
- Ultra-thin

TECHNICAL DATA

- Small Overall Dressing Size: 1.00 x 2.88in. Dressing Pad Size (diameter): 0.63in.
- Medium Overall Dressing Size: 1.50 x 2.63in. Dressing Pad Size (diameter): 1.00in.
- Large Overall Dressing Size: 2.00 x 3.25in. Dressing Pad Size (diameter): 1.50in.
- Acrylic adhesive
- Polyethylene and Polyurethane films

ORDERING INFO

Part Number	Description	Case Size
C013080-SOBC	6-Pack of Bandages - Small	70 pkgs/case
C013080-MOBC	6-Pack of Bandages - Medium	70 pkgs/case
C013080-LOBC	6-Pack of Bandages - Large	70 pkgs/case
C013080-MIXOBC	6-Pack of Bandages - Variety 2 of each size	70 pkgs/case
C013080-100S	100 Count Box - Small	4 boxes/case
C013080-100M	100 Count Box - Medium	4 boxes/case
C013080-100L	100 Count Box - Large	4 boxes/case



YOUR AUTHORIZED SAM MEDICAL PRODUCTS DEALER