

SAM SPLINT MARKETING TEXT

Tagline

The gold standard in splints

Or

For immobilizing skeletal and soft tissue injuries. The gold standard in splints.

Short Description (21 words)

The SAM® Splint is a lightweight, flexible splint that can be bent for strength to splint any fractured or injured limb.

Medium Description (89 words)

The SAM® Splint is the gold standard in splints. Built from a thin core of aluminum alloy and sandwiched between two layers of closed-cell foam, the SAM Splint can be bent into any of three simple curves, becoming extremely strong and supportive for any fractured or injured limb. The SAM Splint is incredibly lightweight (only 4.0 oz.), and no other device offers so many variations for the treatment of orthopedic emergencies. Also available is the SAM Splint XL at 5.5" wide, for improved support and immobilization of larger individuals.

Long Description (145 words)

The SAM® Splint is the gold standard in splints. Built from a thin core of aluminum alloy and sandwiched between two layers of closed-cell foam, the SAM Splint can be bent into any of three simple curves, becoming extremely strong and supportive for any fractured or injured limb. The SAM Splint is incredibly lightweight (only 4.0 oz.), and no other device offers so many variations for the treatment of orthopedic emergencies. Also available is the SAM Splint XL at 5.5" wide, for improved support and immobilization of larger individuals.

The SAM Splint can be found from mountain trails to military combat zones, from outer space to the ocean depths, in all weather conditions and terrains. The SAM® Splint is the emergency splint of choice for EMS personnel, military medics, athletic trainers, wilderness rescuers, outdoor adventurers, first responders, safety engineers, veterinarians, hospitals, ski patrols and more.

You may also use the descriptive text on the product webpage.