

# SAM SUGAR TONG™

Ease and comfort for the immobilization of acute wrist and forearm fractures




Strong and malleable like the **SAM® SPLINT**, the **SAM SUGAR TONG™** is designed for quick-and-easy application when a patient requires temporary immobilization of acute wrist and forearm fractures. The aluminum core allows it to be molded to fit the patient's wrist and forearm. This is ideal for post casting immobilization and post surgical immobilization.

## FEATURES AND BENEFITS

- Limits flexion, extension and rotation of the wrist
- Rapidly applied in less than one minute without water or gloves
- Soft padding provides excellent ventilation and comfort
- Easily re-molded
- Can be easily removed and reapplied
- Lightweight
- Radiolucent
- Washable

## TECHNICAL DATA

- Adult: 4.25" × 37"; 6.2 oz.
- Pediatric: 4.25" × 23"; 4 oz.



CE Approved  
ISO 13485:2003

## ORDERING INFO

Part Number	Description	Case Size
	Black	
ST052106	SAM Sugar Tong™ Splint, Adult 37"	4/case
ST052106P	SAM Sugar Tong™ Splint, Ped. 23"	4/case



**SAM** MEDICAL  
PRODUCTS®

customerservice@sammedical.com  
P.O. Box 3270 Tualatin, OR 97062  
800.818.4726 | 503.639.5474  
www.sammedical.com

