



SAM
MEDICAL

SAM IO TRAINING KITS

TRAIN. PREPARE. PRESERVE LIFE.

SAM IO TRAINING KITS: A VERSATILE WAY OF TRAINING

The SAM IO Training Kits include all components needed to simulate intraosseous access with the SAM IO system. Every kit includes a nylon Training Case, Driver, three Needle Sets, Stabilizer, and Adaptor. Simulation Sawbones, including the proximal humerus, proximal tibia, and distal tibia, are also available in the SAM IO Training Kits. All included parts are for training purposes only.

ENGINEERED FOR SURVIVAL

MANUALLY OPERATED DRIVER

SAM IO is a manually operated intraosseous access system. Catheter placement is achieved by continuously actuating (repeatedly compressing) the driver's trigger assembly while gently guiding the needle assembly into position. Repeated, full trigger actuation creates rotational spin of the needle assembly which, when combined with gentle downward pressure, results in controlled IO placement. Once needle assembly is properly positioned, the stylet is removed to expose standard Luer-lock for extension set connection. With the extension set connected, aspiration verification, flushing and selected treatment(s) may commence.





SAM[®] IO

KITS AVAILABLE:

- 1 PN: IO720-EN**
 x1 Training Case
 x1 Driver
 x1 Stabilizer
 x1 Adaptor
 x1 15 mm Needle Set
 x1 25 mm Needle Set
 x1 45 mm Needle Set
 x1 Proximal Humerus Sawbone

- 2 PN: IO721-EN**
 x1 Training Case
 x1 Driver
 x1 Stabilizer
 x1 Adaptor
 x1 15 mm Needle Set
 x1 25 mm Needle Set
 x1 45 mm Needle Set
 x1 Proximal Tibia Sawbone

- 3 PN: IO722-EN**
 x1 Training Case
 x1 Driver
 x1 Stabilizer
 x1 Adaptor
 x1 15 mm Needle Set
 x1 25 mm Needle Set
 x1 45 mm Needle Set
 x1 Distal Tibia Sawbone

- 4 PN: IO723-EN**
 x1 Training Case
 x1 Driver
 x1 Stabilizer
 x1 Adaptor
 x1 15 mm Needle Set
 x1 25 mm Needle Set
 x1 45 mm Needle Set
 x1 Proximal Humerus Sawbone
 x1 Proximal Tibia Sawbone
 x1 Distal Tibia Sawbone



TRAINING REPLACEMENT NEEDLES

For training & demonstration purposes



PRODUCT SPECIFICATIONS

CASE SIZE & WEIGHT*

Length: 11 in (28 cm)
Height: 5.25 in (13 cm)
Width: 8.5 in (22 cm)
Weight**: 34 oz (963 g)

AVAILABLE KITS

SAM IO Training Kit - IO720-EN
SAM IO Training Kit - IO721-EN
SAM IO Training Kit - IO722-EN
SAM IO Training Kit - IO723-EN

MODELS AVAILABLE


Driver:
IO700-EN
Stabilizer:
1-Pack IO740-1P-EN
5-Pack IO740-5P-EN
Adaptor:
IO742-EN

Needles:

15 mm 1-Pack IO705-1P-EN
15 mm 5-Pack IO705-5P-EN
25 mm 1-Pack IO706-1P-EN
25 mm 5-Pack IO706-5P-EN
45 mm 1-Pack IO707-1P-EN
45 mm 5-Pack IO707-5P-EN

GLOBAL SALES

U.S. Toll Free: 800.818.4726
U.S. Local: 503.639.5474
U.S. Fax: 503.639.5425
sammedical.com

REF IO720-EN	REF IO722-EN
REF IO721-EN	REF IO723-EN
 SAM Medical Products® 12200 SW Tualatin Rd., Ste 200 Tualatin, OR 97062 USA	

* Size and weight vary by configuration and manufacturing process.
** weight represents PN IO720-EN
SAM® is a registered trademark in the United States and other countries.
©2020 SAM Medical Products. All rights reserved.

Training Replacement

Needles:

15 mm - IO725
25 mm - IO726
45 mm - IO727
Mixed Pack - IO728